

Pressure Relieving Gel Supports

Redefining pressure care standards
in the operating theatre



AneticAid

Innovative **Medical** Technology - Practically Applied

Anetic Aid is a world-leading designer and UK-based manufacturer of medical devices, used extensively for patient transport, treatment and recovery.

Our mission

To deliver best outcomes for patients and best experiences for clinical teams: we listen carefully and understand needs. This allows us to design, develop and manufacture – and to establish partnerships – that lead to the supply of innovative, intuitive, high-quality medical devices for patient transport, treatment and recovery.

Pressure Relieving Mattresses and Gel Supports in this brochure reflect these values. In the UK, dealing with pressure sores accounts for about 4% of national healthcare expenditure and the risk of pressure sores during surgery is well documented.

It is also widely recognised that pressure ulcers – and the risk of infection they bring – are often preventable, so having the appropriate equipment and techniques for dealing with surgical patients is a top priority. For over 25 years, Anetic Aid has been the sole UK distributor for Trulife's Oasis, Oasis Elite, Oasis+ and Azure ranges of high quality gel pads and positioning devices..



Contents

K8 Pressure Care Mattresses	3
K8 Technology, Operating Table and Trolley Mattresses	4-5
Gel Technology, Head and Neck	6-7
Arm and Hand	8-9
Leg, Knee and Foot	10-11
Shoulder, Spine and Torso	12-13
Gynaecology	14-15
General	16-17
Patient Positions	18
Purchase Guide	19

 Made to order



PRESSURE CARE MATTRESSES






Eschmann Type

Description	Code
Eschmann T10 / T20 Type K8 Mattress - Complete with Touch Fastener	 10081
Eschmann T10 / T20 / T30 / T50 Type K8 Mattress - Head Section Only with Touch Fastener - 330mm	 10081-A1
Eschmann T10 / T20 Type K8 Mattress - Upper & Lower Torso Section with Touch Fastener - 1110mm	 10081-B1
Eschmann T10 / T20 / T30 / T50 Type K8 Mattress - Leg Section Only with Touch Fastener - 540mm	 10081-C1

Maquet Type

Description	Code
Maquet Alpha / Beta Type K8 Mattress - Head - Tapered 510 to 460mm	 10079-A1
Maquet Alpha / Beta Type K8 Mattress - Shoulder / Foot Section - 250mm	 10079-B1
Maquet Alpha / Beta Type K8 Mattress - Body Section - 860mm	 10079-C1
Maquet Alpha / Beta Type K8 Mattress - Leg Section - 530mm	 10079-D1

Description	Code
Maquet Alpha / Beta Type K8 Mattress - Head Section 300mm	 10079-A2
Maquet Alpha / Beta Type K8 Mattress - Body Section - 990mm	 10079-C3
Maquet Alpha / Beta Type K8 Mattress - Leg - 570mm	 10079-D2



Combining industry knowledge with innovative technology



Specialist equipment constructed using the latest production techniques



Manufacturer, designer and supplier of high quality medical equipment around the world




Science of Kinetics: human body's eight points of contact on a hard surface


K8 Technology

Manufactured from temperature-sensitive gel feel foam that moulds to the patient's body shape, relieving pressure by distributing the load as evenly as possible

Outer material is anti-microbial, water resistant, vapour permeable and latex free for optimum infection control.



K8 OPERATING TABLE AND TROLLEY MATTRESSES



Maquet Orthostar Type

Description	Code
Maquet Orthostar Type K8 Mattress - Head - Tapered 510 to 460mm	✂ 10079-A1
Maquet Orthostar Type K8 Mattress - In-fill - 260mm	✂ 10078-B1
Maquet Orthostar Type K8 Mattress - Upper Torso - 600mm	✂ 10078-C1
Maquet Orthostar Type K8 Mattress - Lower Torso - 440mm	✂ 10078-D1
Maquet Orthostar Type K8 Mattress - Hip In-fill - Left	✂ 10078-E1
Maquet Orthostar Type K8 Mattress - Hip In-fill - Right	✂ 10078-F1
Maquet Orthostar Type K8 Mattress - Leg Section 680mm - PAIR	✂ 10079-E1



Steris Type

Description	Code
Steris Type Operation Table - K8 Mattress	✂ 10086

Merivaara Promerix Type

Description	Code
Merivaara Promerix Type Operation Table - K8 Mattress	✂ 10041



QA4 Surgery Trolley Type

Description	Code
QA4 Surgery Trolley - K8 Mattress - Body Section Only	21360B
QA4 Surgery Trolley - K8 Mattress - Head Section Only	21360H
QA4 Surgery Trolley - K8 Mattress - Leg Section Only	21360L



QA3 Patient Trolley Type

Description	Code
QA3™ Mattress with Waterproof Zip - for QA3™ v2.0 Trolley	21184
QA3™ Replacement Mattress Cover with Waterproof Zip - v2.0	21184-B1

QA3™ Mattress with Waterproof Zip - for QA3™ v3.0 / v4.0 Trolley	21186
QA3™ Replacement Mattress Cover - for QA3™ v3.0 / v4.0	21186-B1

QA3™ Mattress with Waterproof Zip - for QA3™ <800> Trolley	21187
QA3™ Replacement Mattress Cover - for QA3™ <800>	21187-B1



Patient Trolley Type

Description	Code
Stryker Prime / Transport Trolley Type K8 Pressure Care Mattress - with Zip Cover - 660 (26") Wide Patient Platform - 100mm (4") Deep	✂ 10560-B1
Stryker Prime / Transport Trolley Type K8 Pressure Care Mattress - with Zip Cover - 760 (30") Wide Patient Platform - 100mm (4") Deep	✂ 10565-B1



Arjo Huntleigh Lifeguard - 20 / 50 / 55 - Trolley Type K8 Pressure Care Mattress - 90mm Deep	✂ 10580
--	---------



Mattress measuring service available

Gel technology to minimise the risk of pressure sores and infection


Pressure, the amount of force exerted on a given area is often measured in millimetres of mercury (mmHg).

When a force or pressure greater than normal capillary pressure is exerted on a body over time, this can restrict blood flow to the area and cause serious tissue damage.


It is therefore important to reduce this pressure, particularly with the more susceptible parts of the body such as the head in periods of immobility, such as during surgery.

Our gel pads help relieve pressure by increasing the area of contact between the body and the supporting surface.



- Latex free
- Reusable





HEAD AND NECK





Head Protector	
Description	Code
Azure Head Protector (AZ100) - Infant / Adolescent / Adult - 260 x 250 x 50mm	10885



Prone Head Rest	
Description	Code
Oasis Elite Prone Head Rest (EL141) - Adult - 280 x 240 x 140mm	10871
Oasis Elite Prone Head Rest (EL142) - Adolescent - 273 x 220 x 135mm	10872



Supine Head Rest	
Description	Code
Oasis Elite Supine Head Rest (EL151) - Adult - 275 x 237 x 90mm	10873
Oasis Elite Supine Head Rest (EL152) - Adolescent - 275 x 230 x 65mm	10874



Face Pad	
Description	Code
Oasis Face Pad (OA140) - 165 x 80 x 8mm	10849



Head and Neck Support	
Description	Code
Oasis+ Head & Neck Support (OP151) - Adult - 264 x 220 x 64mm	10900


Closed Head Rings	
Description	Code
Oasis Closed Head Ring (OA010) - Neonate / Newborn - Ø90 x 20mm	10801
Oasis Closed Headring (OA011) - Infant / Adolescent - Ø140 x 35mm	10802
Oasis Closed Headring (OA012) - Adolescent / Adult - Ø200 x 45mm	10803

Adult Closed Head Ring	
Description	Code
Oasis Elite Closed Head Ring (EL011) - Infant / Adolescent - Ø140 x 34mm	10861
Oasis Elite Closed Head Ring (EL012) - Adult - Ø200 x 50mm	10862

Open Head Ring	
Description	Code
Oasis Elite Open Head Ring (EL021) - Infant / Adolescent - Ø140 x 34mm	10863
Oasis Elite Open Head Ring (EL022) - Adult - Ø200 x 50mm	10864



Secure positioning for the head and neck




Armboards and Protectors

Flat armboard pads have a slight indentation to cradle the arm and provide effective pressure relief during various surgical procedures.




- Medical grade silicone
- X-ray translucent




Ulnar / Brachial Nerve Protector

Description	Code
Oasis Ulnar / Brachial Nerve Protector (OA061) - Adolescent / Adult - 460 x 210 x 8mm	10829



Arm Retainer Cover

Description	Code
Oasis Arm Retainer Pad (OA100) - 170 x 110 x 18mm	10834



Lateral Park Bench Positioner

Description	Code
Oasis Elite Lateral Park Bench Positioner (EL913)	10995





Oasis Extremity Belt

Description	Code
Oasis Extremity Belt (OA113) - Pair - 315 x 45 x 10mm	10843



Arm Protector

Description	Code
Oasis Elite Contoured Armboard / Knee Pad (EL401) - Adult - 510 x 155 x 40mm	10869



Armboard Pads

Description	Code
Oasis Armboard Pad (OA070) - Infant / Adolescent / Adult - 395 x 115 x 10mm	10830
Oasis Armboard Pad (OA071) - Infant / Adolescent / Adult - 600 x 200 x 10mm	10831
Oasis Armboard Pad (OA072) - Adolescent / Adult - 600 x 125 x 10mm	10832



Arm Protector

Description	Code
Azure Arm Protector (AZ201) - Adolescent / Adult - 460 x 180 x 35mm	10890




Protection for the wrist and elbow, particularly the ulnar nerve

Legs and Heels

The anatomy of the heel – its combination of nerves and lack of soft tissue over the bone – make it one of the most vulnerable parts of the body during surgery, particularly when the patient is in a supine position.



- Easily Cleaned
- High load bearing capacity



LEG, KNEE AND FOOT



Heel Cup	Code
Description Oasis Heel Cup (OA131) - Neonate / Infant / Adolescent / Adult - Pair - 105 x 90 x 43mm	10846

Contoured Arm / Knee Pad	Code
Description Oasis Contoured Armboard / Knee Pad (OA101) - 510 x 150 x 35mm	10833
Oasis Elite Contoured Armboard / Knee Pad (EL402) - Adolescent / Adult - 395 x 140 x 40mm	10870

Dome Positioner	Code
Description Oasis Elite Dome Positioner (EL210) - Adolescent / Adult - 350 x 150 x 75mm	10876

Lateral Leg Positioner	Code
Description Oasis Elite Lateral Leg Positioner (EL800) - Adult - 650 x 400 x 230mm	10992

Heel Pad	Code
Description Oasis Elite Heel Pad (EL130) - Adolescent / Adult - Pair - 180 x 100 x 70mm	10875

Heel Protector	Code
Description Azure Heel Protector (AZ401) - Infant / Adolescent / Adult - 190 x 175 x 30mm	10891

Leg Protector	Code
Description Azure Leg Protector (AZ400) - Adolescent / Adult - 515 x 180 x 40mm	10892

Heel and Ankle Protector	Code
Description Oasis Heel & Ankle Protector (OA200) - Neonate / Infant / Adolescent / Adult - 1250 x 120 x 10mm	10853

Protective pressure care for vulnerable extremities like the heel



SHOULDER, SPINE AND TORSO

Trulife Oasis Elite option: light weight and pre-contoured

One of four pressure care construction options from Trulife, the Oasis Elite range combines silicone gel and light weight foam with a pre-contoured surface to provide superior pressure relief and patient protection in theatre.

Anatomical pre-contouring facilitates greater conformity between gel and body shape; weight is then redistributed over a larger area.

Oasis Elite products are more than 50% lighter than comparable gel-only products.




- Foam and silicone gel combination
- Enhanced protection
- Light weight
- Highly durable
- Easy to clean and maintain



Montreal Laminectomy Support

Description	Code
Elite Foam Large Prone Positioner (EL904) - Montreal Laminectomy - Adult - 1115 x 520 x 140mm	☞ 10990




Prone Chest Positioner

Description	Code
Oasis Elite Prone Positioner (EL903) - Adult - 610 x 530 x 130mm	☞ 10883




Dome Positioner

Description	Code
Oasis Elite Prone Positioner (EL905)	☞ 10884




Elite Prone Positioning System

Description	Code
Elite Foam ARDS Prone Positioning System (ELP200)	☞ 10991
Elite Foam Disposable Prone Head Positioner (ELP141) Supplied in packs of 3	☞ 10991-A1





Oasis Chest Rolls

Description	Code
Oasis Flat Bottom Chest Roll (OA211) - Adult - 305 x 95 x 75mm	10854
Oasis Flat Bottom Chest Roll (OA212) - Adult - 405 x 95 x 75mm	10855
Oasis Flat Bottom Chest Roll (OA213) - Adult - 505 x 92 x 75mm	10856



Elite Chest Rolls

Description	Code
Oasis Elite Flat Bottom Chest Roll (EL214) - Adult - 300 x 140 x 110mm	☞ 10877
Oasis Elite Flat Bottom Chest Roll (EL215) - Adult - 400 x 140 x 110mm	☞ 10878
Oasis Elite Flat Bottom Chest Roll (EL216) - Adult - 505 x 140 x 110mm	☞ 10879
Oasis Elite Flat Bottom Chest Roll (EL217) - Adult - 405 x 150 x 150mm	☞ 10880
Oasis Elite Flat Bottom Chest Roll (EL218) - Adult - 510 x 150 x 150mm	10881



Patient Positioner

Description	Code
Azure Universal Positioner (AZ500) - Adolescent / Adult - 425 x 155 x 60mm	☞ 10893



Secure supine, lateral or prone positioning

Trulife Azure option:
redefining pressure care


Azure's revolutionary gel combination works together to replicate and act as an additional layer of protective tissue and skin.

It conforms to the individual's shape in a more contoured and less structured manner. This results in effective protection against nerve damage as well as pressure sore prevention.

Azure: advanced construction for state-of-the-art pressure care




- Revolutionary new gel combination
- Contoured protection
- Redefining pressure care
- Advanced surgical positioning and protection
- 12 month guarantee




Crutch Stirrup Pads

Description	Code
Oasis Crutch Stirrup Pad (OA090) - Adolescent / Adult - Pair - 320 x 220 x 10mm	10838
Oasis Crutch Stirrup Pad (OA091) - Adolescent / Adult - Pair - 395 x 245 x 10mm	10839

Crutch Stirrup Pad

Description	Code
Oasis Crutch Stirrup Pad (OA092) - Adolescent / Adult - Pair - 520 x 300 x 10mm	10840







Table Pads

Description	Code
Oasis 3-piece Operating Table Pad (OA037) - Torso Section with Cut-out - 900 x 500 x 10mm	10813
Oasis Perineal Table Pad (OA040) - 520 x 520 x 10mm	10816





Profile Pad – Universal

Description	Code
Azure Profile Pad (AZ610) - Universal - 520 x 520 x 10mm	10886


Profile Pad – ¾ Length

Description	Code
Azure Profile Pad (AZ611) - ¾ Length - 1150 x 520 x 10mm	10887



Sacral Protector

Description	Code
Azure Sacral Protector (AZ300) - Adolescent / Adult - 400 x 380 x 25mm	10888



Versatile pressure care pads for a variety of procedures



Trulife Oasis option:
the original silicone solution

One of four pressure care construction options from Trulife, the Oasis range consists of reusable 100% silicone gel pads.

The silicone is softer than both skin and underlying tissue which means that pressure is relieved by the slight movement of the silicone gel and resulting dissipation of pressure/force across the product.

Oasis: gel pressure care for positioning across a wide range of surgical specialisms



- Silicone superior pressure relief
- Conforms to body shape
- Comprehensive range
- Can be heated or cooled
- Latex free and skin compatible





Patient Positioning

Description	Code
Oasis Universal Patient Positioner (OA050) - Infant / Adolescent / Adult - 150 x 60 x 40mm	10823
Oasis Universal Patient Positioner (OA051) - Infant / Adolescent / Adult - 330 x 70 x 50mm	10824
Oasis Universal Patient Positioner (OA052) - Adolescent / Adult - 290 x 100 x 60mm	10825
Oasis Universal Patient Positioner (OA053) - Adolescent / Adult - 345 x 120 x 70mm	10826




Oasis Patient Positioning Belt

Description	Code
Oasis Patient Positioning Belt (OA110) - 505 x 85 x 8mm	10842




Oasis Non-slip Discs

Description	Code
Oasis+ Non-slip Disc (OP015) - Ø200 x 4mm	10975



Operating Table Pads

Description	Code
Oasis Standard Operating Table Pad (OA030) - Universal - 520 x 520 x 10mm	10807



Operating Table Pads

Description	Code
Full Body - Oasis Standard Operating Table Pad (OA031) - ¾ Length - 1150 x 520 x 10mm	10808
Oasis Standard Operating Table Pad (OA032) - Full Length - 1800 x 520 x 10mm	10809



Vacuum Beanbag Positioners

Description	Code
Vacuum Beanbag Positioner - Full Body - 1000 x 900mm	10275
Vacuum Beanbag Positioner - Square - 500 x 500mm	10276
Vacuum Beanbag Positioner - Head Stabiliser - 400 x 300mm	10278
Vacuum Beanbag Positioner - Full Length - 1650 x 720mm	10279



Exclusive distribution of Trulife products across the UK

SUPINE position

- Administration of General Anaesthesia
- Abdominal surgery
- Open heart surgery
- Surgery of the face, neck or mouth
- Most surgery on extremities



The patient is on their back with their face up towards the ceiling. Their legs are uncrossed and their arms are at their sides, or on arm boards.

(Includes Trendelenburg and reverse Trendelenburg)

PRONE position

- Surgeries to the back, spine and the backs of the legs



The patient is on their abdomen and their arms are at their sides or on armboards.

(Includes Kraske 'Jack Knife' position)

LATERAL position

- Surgery of the chest or lungs
- Surgery on the kidneys
- Hip surgery



The patient is on their side, with their bottom leg flexed and their top leg straight.

The bottom arm is supported on an arm board, and the top arm has special support.

LITHOTOMY position

- Gynaecological surgery
- Genitourinary (GU) surgery



The patient is on their back, their legs are lifted in stirrups and their arms are usually on arm boards.

Code Description

Code	Description	SUPINE	PRONE	LATERAL	LITHOTOMY
10802	Oasis Closed Heading (OA011) - Infant / Adolescent - Ø140 x 35mm	✓			✓
10803	Oasis Closed Heading (OA012) - Adolescent / Adult - Ø200 x 45mm	✓			✓
10807	Oasis Standard Operating Table Pad (OA030) - Universal - 520 x 520 x 10mm	✓	✓	✓	✓
10808	Oasis Standard Operating Table Pad (OA031) - ¾ Length - 1150 x 520 x 10mm	✓	✓	✓	✓
10809	Oasis Standard Operating Table Pad (OA032) - Full Length - 1800 x 520 x 10mm	✓			✓
10813	Oasis 3-piece Operating Table Pad (OA037) - Torso Section with Cut-out - 900 x 500 x 10mm	✓			✓
10824	Oasis Universal Patient Positioner (OA051) - Infant / Adolescent / Adult - 330 x 70 x 50mm	✓	✓	✓	
10825	Oasis Universal Patient Positioner (OA052) - Adolescent / Adult - 290 x 100 x 60mm	✓	✓	✓	
10826	Oasis Universal Patient Positioner (OA053) - Adolescent / Adult - 345 x 120 x 70mm	✓	✓	✓	
10829	Oasis Ulnar / Brachial Nerve Protector (OA061) - Adolescent / Adult - 460 x 210 x 8mm	✓	✓	✓	✓
10831	Oasis Armboard Pad (OA071) - Infant / Adolescent / Adult - 600 x 200 x 10mm	✓	✓	✓	✓
10838	Oasis Crutch Stirrup Pad (OA090) - Adolescent / Adult - Pair - 320 x 220 x 10mm	✓	✓	✓	✓
10839	Oasis Crutch Stirrup Pad (OA091) - Adolescent / Adult - Pair - 395 x 245 x 10mm	✓	✓	✓	✓
10853	Oasis Heel & Ankle Protector (OA200) - Neonate / Infant / Adolescent / Adult - 1250 x 120 x 10mm	✓		✓	
10856	Oasis Flat Bottom Chest Roll (OA213) - Adult - 505 x 92 x 75mm	✓	✓	✓	
10861	Oasis Elite Closed Head Ring (EL011) - Infant / Adolescent - Ø140 x 34mm	✓			✓
10862	Oasis Elite Closed Head Ring (EL012) - Adult - Ø200 x 50mm	✓			✓
10863	Oasis Elite Open Head Ring (EL021) - Infant / Adolescent - Ø140 x 34mm		✓	✓	
10864	Oasis Elite Open Head Ring (EL022) - Adult - Ø200 x 50mm		✓	✓	
10869	Oasis Elite Contoured Armboard / Knee Pad (EL401) - Adult - 510 x 155 x 40mm	✓	✓	✓	✓
10871	Oasis Elite Prone Head Rest (EL141) - Adult - 280 x 240 x 140mm		✓		
10873	Oasis Elite Supine Head Rest (EL151) - Adult - 275 x 237 x 90mm	✓			
10875	Oasis Elite Heel Pad (EL130) - Adolescent / Adult - Pair - 180 x 100 x 70mm	✓			
10876	Oasis Elite Dome Positioner (EL210) - Adolescent / Adult - 350 x 150 x 75mm		✓	✓	
10881	Oasis Elite Flat Bottom Chest Roll (EL218) - Adult - 510 x 150 x 150mm	✓	✓	✓	
10882	Oasis Elite Lateral Positioner (EL900) - Adolescent / Adult - 720 x 500 x 160mm		✓		
10884	Oasis Elite Prone Positioner (EL905) - 675 x 500 x 100mm		✓		
10885	Azure Head Protector (AZ100) - Infant / Adolescent / Adult - 260 x 250 x 50mm	✓			✓
10886	Azure Profile Pad (AZ610) - Universal - 520 x 520 x 10mm	✓	✓	✓	✓
10887	Azure Profile Pad (AZ611) - ¾ Length - 1150 x 520 x 10mm	✓	✓	✓	✓
10888	Azure Sacral Protector (AZ300) - Adolescent / Adult - 400 x 380 x 25mm				✓
10890	Azure Arm Protector (AZ201) - Adolescent / Adult - 460 x 180 x 35mm	✓	✓	✓	✓
10891	Azure Heel Protector (AZ401) - Infant / Adolescent / Adult - 190 x 175 x 30mm	✓			
10892	Azure Leg Protector (AZ400) - Adolescent / Adult - 515 x 180 x 40mm	✓			
10990	Elite Foam Large Prone Positioner (EL904) - Montreal Laminectomy - Adult - 1115 x 520 x 140mm		✓		
10992	Oasis Elite Lateral Leg Positioner (EL800) - Adult - 650 x 400 x 230mm			✓	



Scan for Trulife pressure care training video



DESIGN & INNOVATION TECHNICAL EXCELLENCE GLOBAL DISTRIBUTION

OUR PRODUCT AREAS

Operating Table Accessories

Patient positioning solutions for you and your patients



SEE ALSO

Operating Table Accessories:

Patient positioning solutions for you and your patients

PLUS
our other product areas



QA3
Patient Trolley System

Compact, light weight and manoeuvrable – the complete solution for patient transport

QA3
Emergency Trolley

Radio-transparent platform with exceptional X-ray access and low height capability

QA4
Surgery Trolley System

The complete multi-purpose mobile operating system

Tourniquet System

Ideal for single or bi-lateral orthopaedic surgery and pain management

Stainless Steel Theatre Furniture

Quality materials and workmanship guaranteed

AnetiCare Service & Maintenance

Product knowledge combined with customer care

For more information on any of our products or service contracts, call:
T +44 (0)1943 878647 sales@aneticaid.com aneticaid.com



Innovative medical technology – practically applied



Anetic Aid Ltd, Sapper Jordan Rossi Business Park, Otley Road, Baildon BD17 7AX, UK

