

SEE
ALSO



QA3 Ophthalmic Stretcher

Outstanding surgical access –
exceptional patient comfort

PLUS

Our other product areas



OUR PRODUCT AREAS



Scan this code to
start a **WhatsApp** chat
with Anetic Aid

Operating Theatre Accessories

Innovative solutions for you and your patients



For more information on any of our products or service contracts, call:
sales@aneticaid.com **+44 (0)1943 878647** aneticaid.com



Innovative **Medical** Technology - **Practically Applied**



Anetic Aid Ltd, Sapper Jordan Rossi Business Park, Otley Road, Baildon BD17 7AX, UK

Innovative **Medical** Technology - **Practically Applied**



Innovative **Medical** Technology -
Practically Applied

HEAD AND NECK

Anetic Aid is a world-leading designer and UK-based manufacturer of medical devices, used extensively for patient transport, treatment and recovery.

Our mission

To deliver best outcomes for patients and best experiences for clinical teams: we listen carefully and understand needs.

This allows us to design, develop and manufacture – and to establish partnerships – that lead to the supply of innovative, intuitive, high-quality medical devices for patient transport, treatment and recovery.

The wide range of patient positioning and pressure care products in this brochure reflect these values. Designed and produced to the highest quality, they can be used with any operating table or surgical stretcher – as well as our own QA4 Mobile Surgery System.

The focus is on providing patient comfort and safety combined with the practical requirements of diverse surgical procedures. The same applies to our Stainless Steel Furniture range – functionality combined with durability. Some specialist items are made to order, while many others are held in stock, for rapid delivery around the country.

Contents

| PATIENT POSITIONING | |
|---------------------------|-------|
| Head and Neck | 3 |
| Arm and Hand | 4-6 |
| Leg, Knee and Foot | 7 |
| Shoulder, Spine and Torso | 8-10 |
| Gynaecology | 11-13 |
| General | 14-17 |

| PRESSURE RELIEVING GEL | |
|---------------------------|-------|
| Head and Neck | 18-19 |
| Arm and Hand | 20-21 |
| Leg, Knee and Foot | 22-23 |
| Shoulder, Spine and Torso | 24-25 |
| Gynaecology | 26-27 |
| General | 28-29 |
| Patient Positions | 30 |
| Purchase Guide | 31 |

| STAINLESS STEEL FURNITURE | |
|--|----|
| Instrument Trolleys | 32 |
| Dressing / Equipment Trolleys | 33 |
| Surgeon's Seating | 34 |
| Transfusion / Infusion Stands – Mobile | 35 |
| Mayo Tables and Bowl Stands | 36 |
| Kick-about and Surgeon's Platform | 37 |
| General | 38 |
| QA3 & QA4 Stretchers | 39 |

✂ Made to order ⚙ Factory fit – please specify at time of ordering



Rubens
Pillow

Ideal for ophthalmics, providing secure positioning of the head while allowing the surgeon to rest their elbows on the pillow's 'wings'.

| Description | Code |
|----------------------|-------|
| Rubens Pillow, Adult | 10180 |



Headrests

Both headrests offer K8 Pressure Care, and the Whitlock Horseshoe also provides special weighting and adjustability.



| Description | Code |
|-----------------------------|-------|
| Doughnut Head Ring, Adult | 10170 |
| Whitlock Horseshoe Headrest | 10185 |



Vacuum
Beanbag

One of our full range of stabilising vacuum beanbags, this model is specially shaped for head positioning.

| Description | Code |
|---|-------|
| Vacuum Beanbag Positioner, Head Stabiliser 400 x 300mm | 10278 |




Theatre
Pillow

Resistant to fluids, but breathable for patient comfort, the pillow has a specialist gas permeable valve to guard against accidental movement when air is expelled.

| Description | Code |
|----------------|-------|
| Theatre Pillow | 10545 |


ARM AND HAND



10110

Arm Retainers and Retainer Pads

| Description | Code |
|--|-------|
| Arm Retainer, L Shaped 235 x 235 x 127mm | 10105 |
| Arm Retainer, J Shaped 235 x 210 x 127mm | 10110 |
| Arm Retainer, L Shaped 235 x 235 x 300mm | 10112 |
| Arm Retainer Pad | 10115 |
| Arm Retainer 300mm Wide Pad | 10117 |



10105 / 10115



10220

Narrow Arm Tables


A traditional and affordable arm table for a wide range of procedures.

| Description | Code |
|---------------------------------------|--------|
| Narrow Arm Table 590 x 210mm c/w Pad | *10220 |
| Narrow Arm Table 590 x 210mm Pad ONLY | 10221 |

*Requires 1 x Clamp (Code: 10250)



10221




10774


Easy Armboard

Attaches via a direct-on mechanism and easy to operate.

| Description | Code |
|--|-------|
| Easy Armboard Lift-to-Position for EU Side Bar, 75mm c/w Pad | 10774 |
| Easy Armboard Lift-to-Position for UK Side Bar, 75mm c/w Pad | 10775 |
| Easy Armboard Lift-to-Position Replacement Pad, 75mm | 10705 |



10775




10843


Extremity Belt and Wristlet with Strap

Useful for many procedures, this product provides extra security when maintaining the position of the extremities on the operating table.

| Description | Code |
|---|------------|
| Oasis Extremity Belt (OA113) - Pair - 315 x 45 x 10mm | 10843 |
| Wristlet with Strap | ✂ 10460-A2 |




10843



10425

Reison Arm and Hand Surgery Table


| Description | Code |
|---|----------|
| Arm & Hand Surgery Table - Carbon Fibre (10-396) - Hourglass - with Twin Hook-on Clamps - includes Pad - L820 x W525 x H720 to 1150mm - for UK Side Bar | 10425-UK |
| Arm & Hand Surgery Table - Carbon Fibre (10-396) - Hourglass - with Twin Hook-on Clamps - includes Pad - L820 x W525 x H720 to 1150mm - for EU Side Bar | 10425-EU |



10425

ARM AND HAND





10233

Arm and Hand Surgery Tables

Anetic Aid offers a wide selection of arm supports including arm tables for positioning and protecting patients in prone, supine or lateral positions.

One of these is our extended ultra light weight arm table which attaches directly, without the need for clamps.

The telescopic leg provides additional stability – ideal for delicate hand surgery or where surgeons or practitioners will need to put weight on the arm table during the procedure.

| Description | Code |
|---|---------|
| Extended Arm Table, X-ray Translucent c/w Pad | 10233 |
| Extended Arm Table 75mm Pad ONLY | ✂ 10234 |




- Extended Arm Table:**
- Attaches directly without the need for clamps
 - X-ray translucent
 - Removable telescopic leg

ARM AND HAND

Reison Rotating Armboard

Operated via a highly secure one-handed ball joint release mechanism, this armboard offers accurate positioning, patient comfort and excellent surgical access.



*Requires 1 x Cushion / Pad (Code: 10415-B1)



A tried and tested armboard which is easily manoeuvred into the desired position without stress to either the shoulder or elbow joint.



| Description | Code |
|--|-----------|
| Armboard - Rotatable (10-387) - with Rotating Clamp for UK Side Bar | *10400-UK |
| Armboard - Rotatable (10-387) - with Rotating Clamp for EU Side Bar | *10400-EU |
| Armboard Cushion / Pad (10-382-H) - Medisoft (Viscoelastic) - Tapered Bottom with Guttered Top Surface - 40mm Deep | 10415-B1 |



Carter Braine Arm Support

This Carter Braine Arm Support is a useful and effective accessory, providing comfortable and convenient placement for the operative or non-operative arm.

| Description | Code |
|------------------------------------|--------|
| Carter Braine Arm Support c/w Pad | *10195 |
| Carter Braine Arm Support Pad ONLY | 10196 |

*Requires 1 x Clamp (Code: 10250)




Lateral / Prone Arm Support

This malleable handrest enables comfortable and convenient positioning of the patient's non-operative arm in any desired orientation, to prevent potential injury.

| Description | Code |
|-------------------------------------|--------|
| Lateral / Prone Arm Support c/w Pad | *10773 |

*Requires 1 x Clamp (Code: 10248)

LEG, KNEE AND FOOT



Tibia Support

This Tibia Support is designed to support the knee or the foot of a patient during knee or leg surgery.

| Description | Code |
|---|-----------|
| Tibia Support (10-367) - L270 x Ø135mm with L770mm Stem | *10440-A1 |

*Requires 1 x Clamp (Code: 10250)



Arthroscopic Stress Post

This positioner provides a levering fulcrum, designed for improved access to the medial compartment of the knee during arthroscopy procedures.

| Description | Code |
|---|--------|
| Arthroscopic Stress Post c/w Fitted Pad | *10776 |

*Requires 1 x Clamp (Code: 10248)




Knee Stabiliser

A convenient, stable device to support the flexed leg during total knee procedures, this 300mm x 100mm foot support eliminates awkward sandbags and simply mounts to the table rail.

| Description | Code |
|-------------------------|--------|
| Knee Stabiliser c/w Pad | *10785 |

*Requires 1 x Clamp (Code: 10250)



Abduction Board

K8 Pressure Care extension pads for stretchers, trolleys or operating tables: helps protect a patient's heels after leg abduction for procedures such as Varicose Vein removal.

| Description | Code |
|---|-------|
| Abduction Board (varicose vein) 75mm Mattress c/w Pads | 10155 |
| Abduction Board (varicose vein) 75mm Mattress Pads ONLY | 10156 |

Positioning for the operative or non-operative arm



DESIGN & INNOVATION



TECHNICAL EXCELLENCE



GLOBAL DISTRIBUTION

SHOULDER, SPINE AND TORSO

Lateral Brace 3-Piece

As well as being appropriate for most makes of operating table, these accessories can be used with Anetic Aid's QA4 Day Surgery System – the complete mobile operating system.

| Description | Code |
|--|---------|
| Lateral Brace - 3-piece Set - includes Two Rectangular Supports & One Circular Support - Fitted Pads | **10770 |
| Lateral Brace - Circular Support Only - Fitted Pad | *10771 |
| Lateral Brace - Rectangular Support Only - Fitted Pad | *10772 |
| Lateral Brace - Replacement Pad - Circular Support | 691120 |
| Lateral Brace - Replacement Pad - Rectangular Support - with Centred Post Hole | 691125 |
| Lateral Brace / Arthroscopic Stress Post - Pad - Off-set | 691130 |

**Requires 3 x Clamps (Code: 10761)
*Requires 1 x Clamp (Code: 10761)



Reison Lateral Supports

This type of support includes two pads, one flat and one curved to accommodate different patient physiques and can be fixed in a number of positions.

| Description | Code |
|--|-------------|
| Lateral Support - Narrow - Rotatable (10-560) - includes one straight & one curved pad | ✂ *10430-L1 |
| Lateral Support - Wide - Rotatable (10-561) - includes one straight & one curved pad | ✂ *10430-L2 |

*Requires 1 x Clamp (Code: 10250)



Lateral / Pelvic Supports

The lateral position is where the patient has their non-operative side on the surgical surface. Our Lateral Support is designed specifically for this position and fits into the small of the back to safeguard the patient from rolling off the operating table.

Our cylindrical, upright Pelvic Support is also designed to position the patient securely in the lateral position. Normally used in conjunction with the Lateral Support, this accessory is sited in the patient's groin / pelvic area.

| Description | Code |
|--------------------------|--------|
| Lateral Support c/w Pad | *10200 |
| Lateral Support Pad ONLY | 10201 |
| Pelvic Support c/w Pad | *10205 |
| Pelvic Support Pad ONLY | 10206 |

*Requires 1 x Clamp (Code: 10250)



SHOULDER, SPINE AND TORSO



Montreal Laminectomy Support

The Montreal Laminectomy Support is a positioning device used in spinal procedures – and any other surgeries requiring the prone position.

The central hollow relieves pressure on the diaphragm ensuring that the patient can breathe as normally as possible.

Manufactured to K8 pressure care specification and ergonomically designed to ensure maximum patient comfort, all edges are sealed for easy cleaning.

| Description | Code |
|--------------------------------------|---------|
| Montreal Laminectomy Support | 10160 |
| Montreal Laminectomy Support - Large | ✂ 10162 |



Shoulder Supports

This shoulder support system has been specially designed to secure the patient in the Trendelenburg position.



It supports without direct pressure on the brachial plexus, protecting against nerve damage. With thick foam padding and a flexible cover, the system supports a patient weighing up to 227kg.

| Description | Code |
|---|--------|
| Shoulder Supports c/w Pads | *10787 |
| Shoulder Supports Replacement Pads (left and right) | 10706 |

*Requires 2 x Clamps (Code: 10761)



Positioning for total hip replacement:
accurate alignment for optimum surgical outcomes



Hip Fix™ has more than
20 years of clinical history

Total hip replacement is an exacting surgical technique and it is important for the orthopaedic surgeon to be able to accurately orientate the implant components.

The Hip Fix™ has a three-point fixation for the lateral position, securing both sides of the iliac crest plus the sacrum, giving the surgeon maximum access while maintaining total stability. Although the position is rigidly maintained, there is no external pressure on the abdomen, ensuring greater patient comfort.

| Description | Code |
|--|----------|
| Hip Fix™ (2015-200) 3-point Fixation System - includes sacrum positioner & vertically mounted iliac crest positioners - disposable pads | 10470 |
| Hip Fix™ (1915-200) 3-point Fixation System - individual support adjustment - includes sacrum positioner & individual vertically mounted iliac crest positioners - disposable pads | 10470-A2 |
| Hip Fix™ Rotary Clamp (10-305-P) - Hook-on for UK Side Bar with Flat Pushing Plunge | 10470-C1 |
| Hip Fix™ Rotary Clamp (10-305-P) - Hook-on for EU Side Bar with Flat Pushing Plunge | 10470-C2 |
| Hip Fix™ Single-use Pad Set - includes one sacrum & two iliac crest positioner pads (Pack of 25) | 10470-P1 |

*Requires 2 x Clamps (Code: 10470-C1) or 2 x Clamps (Code: 10470-C2)

Single-use pads

The Hip Fix's stabilising arms are fitted with single-use pads for each new procedure, reducing the risk of cross-contamination – vital for such an invasive surgical procedure.

- Single use only
- Latex-free and anti-allergenic
- Simple, hygienic and safe



Finned Elite Stirrups

The unique design of the Finned Elite model Lift-Assist™ Stirrup, with its 'fin'-shaped support at the knee and squeeze grip handle, offers even greater control when manoeuvring patient weights of up to 227kg.

| Description | Code |
|---|--------|
| Lift-Assist™ Lithotomy Stirrups - Finned Elite - includes Boot-pads | *10757 |

*Requires 2 x Clamps (Code: 10761)

Baxter

Accessory Clamps supplied separately unless otherwise stated

Three-point fixation for the lateral position



DESIGN & INNOVATION



TECHNICAL EXCELLENCE



GLOBAL DISTRIBUTION

GYNAECOLOGY

Lift-Assist™ Lithotomy Stirrups

The lithotomy position is widely used for procedures ranging from simple pelvic examinations to surgeries in gynaecology, urology, and other gastrointestinal specialisms.

The advantages for the medical practitioner are obvious in terms of visibility and access, but care has to be taken to protect the patient from potential nerve damage caused by pressure, and from over extension of abductor muscles.

Protecting patients

Our boot type Lift-Assist™ Lithotomy Stirrups can significantly help to reduce risks for patients: sold in pairs, all of our stirrups have self adjusting boots designed to cover and protect the head of the fibula and hence, the peroneal nerve.

They minimise pressure on the calf when moving the stirrup and the hydraulic mechanism allows easy movement of the leg when placing it in the desired position without compromising the sterile field. Simply releasing the twist grip handle secures the leg holder in all directions.



| Description | Code |
|---|--------|
| Lift-Assist™ Lithotomy Stirrups - Standard - includes Boot-pads | *10755 |
| Lift-Assist™ Lithotomy Stirrups - Heavy Duty - includes Boot-pads | *10758 |

*Requires 2 x Clamps (Code: 10761)

Quick reference guide:

| STANDARD | FINNED ELITE | HEAVY DUTY |
|--|--|--|
| For patient weights up to 159 kg | For patient weights up to 227 kg | For patient weights up to 327 kg |
| Handle type: Twist Grip | Handle type: Squeeze Grip | Handle type: Twist Grip |
| Boot type: Standard | Boot type: Protective Fin | Boot type: Enlarged Standard |
| Lithotomy range: +84° / -33° | Lithotomy range: +84° / -33° | Lithotomy range: +84° / -33° |
| Abuction range: +25° / -9° | Abuction range: +25° / -9° | Abuction range: +25° / -9° |
| Code: 10755 | Code: 10757 | Code: 10758 |

GYNAECOLOGY

Mobile Storage Cart

This is the perfect way to store and protect your investment: a strong and light weight cart which holds and transports stirrups or legholders.

Compact and state-of-the-art, it can take stirrups from any manufacturer as it has two different size hooks (up to a 19mm rod on the smaller hook and up to a 25mm rod on the larger one).

The sturdy, moulded, four-caster polymer base is dent-resistant and the deep storage bin is a great place to store clamps and pads securely.

| Description | Code |
|--|-------|
| Storage Cart For Lift-Assist™ Lithotomy Stirrups | 10760 |

Boot Pad Range

| Description | Code |
|--|--------|
| Boot Pad for Standard and Elite Stirrups - 10755 / F-LH350 + 10756 / F-LH500 | 691200 |
| Boot Pad for Heavy Duty Stirrups - 10758 / F-LH720 | 691205 |
| Boot Pad for PAL® / PAL Pro® Stirrups 10700 / A-10023 + 10750 / A-10052 | 691210 |
| Boot Pad for Yellofins® Stirrups 10752 / O-YFASI + 52938 / O-YFES | 691215 |
| Boot Pad for Ultrafins® Stirrups - 53411 / O-UFASE | 691220 |

Knee Crutch Stirrups

| Description | Code |
|---|--------|
| Simplicity™ Knee Crutch Stirrups c/w Boot Pads (pair) | *10765 |
| Boot Pads - for Simplicity™ Knee Crutch Stirrups - Pair | 10766 |

*Requires 2 x Clamps (Code: 10248)

Cysto-Lift Knee Crutch Stirrups

| Description | Code |
|---|-----------|
| Cysto-Lift Knee Crutch - includes Pads (pair) | *10765-A1 |
| Boot Pads - for Simplicity™ Knee Crutch Stirrups - Pair | 10766 |

*Requires 2 x Clamps (Code: 10761)

Clamps

We offer a full range of clamps, both direct-on and fixed for UK and EU Side Bar.

| Description | Code |
|--|-------|
| Clamp - Basic Blade Handwheel - for UK and EU Side Bar | 10761 |
| Clamp, Clip-on™ Blade Handwheel - for UK Side Bar | 10763 |
| Clamp, Clip-on™ Blade Handwheel - for EU Side Bar | 10764 |



Vacuum Beanbags

Made from strong and pliable polyurethane and filled with thousands of tiny plastic beads, air is extracted to keep the patient securely positioned for the duration of the procedure.

| Description | Code |
|--|-------|
| Vacuum Beanbag Positioner - Full Body - 1000 x 900mm | 10275 |
| Vacuum Beanbag Positioner - Square - 500 x 500mm | 10276 |
| Vacuum Beanbag Positioner - Full Length - 1650 x 720mm | 10279 |




Table Width Extenders

| Description | Code |
|---|---------------|
| Operation Table Widening Section - EU - L240mm - Pair | ✂ 10450-EU-A1 |
| Operation Table Widening Section - EU - L300mm - Pair | ✂ 10450-EU-A2 |
| Operation Table Widening Section - EU - L445mm - Pair | ✂ 10450-EU-A3 |

| | |
|---|------------|
| Operation Table Widening Section - EU - L510mm - Single | ✂ 10788-EU |
| Operation Table Widening Section - UK - L510mm - Single | ✂ 10788 |



Security Strap

| Description | Code |
|---|-------|
| Security Strap, Medium c/w Padded Centre Section NO HOOKS | 10792 |
| Security Strap, Large c/w Padded Centre Section NO HOOKS | 10793 |



Airway Management

| Description | Code |
|---------------------------------|-------|
| Anaesthetic Tube Holder, Flat | 10120 |
| Anaesthetic Tube Holder, Angled | 10122 |
| Flexible Tube Holder - Large | 10124 |



Anaesthetic Screen


A traditional support for anaesthetic screens.

| Description | Code |
|----------------------------|--------|
| Anaesthetic Screen c/w Pad | *10190 |
| Anaesthetic Screen Sleeve | 10192 |

*Requires 1 x Clamp (Code: 10250)

Anaesthesia Screens

An extended version of the anaesthetic screen with capacity for broader drapes.



| Description | Code |
|--|------------|
| Anaesthesia Frame - Carbon Fibre (10-313) - Hook-on - UK | ✂ 10460-A4 |
| Anaesthesia Frame - Carbon Fibre (10-313) - Hook-on - EU | ✂ 10460-A5 |
| Anaesthesia Screen - Flexible | ✂ *10784 |

*Requires 1 x Clamp (Code: 10248)



Oxygen Delivery Bar and Drape Screen

| Description | Code |
|------------------------------------|--------|
| Oxygen Delivery Bar & Drape Screen | *21365 |

*Requires 1 x Clamp (Code: 10250)

Patient security on the operating table



DESIGN & INNOVATION



TECHNICAL EXCELLENCE



GLOBAL DISTRIBUTION



Needle Triever™ Device

It doesn't happen often, but if a suture needle drops onto the floor during a procedure, it's obviously essential to find it – and find it quickly.

Our magnetic 'Needle Triever' is the perfect solution, reaching into awkward corners or under equipment – especially the operating table – easily and effectively.

| Description | Code |
|------------------------|-------|
| Needle Triever™ Device | 10796 |

Sandbags

Positioning stability – outer covers made of the same easy-to-clean, breathable membrane used in our mattress range




| Description | Code |
|---|-------|
| Sandbag - Rectangular - L200 x W150mm SMALL | 10255 |
| Sandbag - Rectangular - L300 x W200mm MEDIUM | 10256 |
| Sandbag - Rectangular - L380 x W250mm LARGE | 10257 |
| Sandbag - Rectangular - L450 x W300mm EXTRA LARGE | 10258 |

D Shapes



| Description | Code |
|---|-------|
| Positioner - D-Shaped (Knee) - W550 x D300 x H175mm | 10295 |

Wedges



| Description | Code |
|---|-------|
| Crawford Wedge Small - W550 x D300 x H100mm | 10550 |
| Crawford Wedge Large - W550 x D510 x H200mm | 10555 |



10248



10250



10761



10763

Clamps

We offer a full range of clamps, both direct-on and fixed for UK and EU Side Bar.

| Description | Code |
|---|-------|
| Clamp - Rotary Socket - Stainless Steel - Hook-on - Universal for UK / EU / US Side Bar | 10248 |
| Universal Clamp, Socket and Rail for UK and EU Side Bar | 10250 |
| Clamp, Basic Blade Handwheel for UK and EU Side Bar | 10761 |
| Clamp, Clip-on™ Blade Handwheel for UK Side Bar | 10763 |
| Clamp, Clip-on™ Blade Handwheel for EU Side Bar | 10764 |

Pre-operative Glove



A safe and hygienic way of removing hair from patients after shaving or clipping, reducing the risk of surgical site contamination and post-procedural infections.

| Description | Code |
|---|-------|
| Pre-operative Glove, Adhesive - Single-use (box of 500) | 11000 |

QA4 Accessory Stand



Accessories need storage when not in use – and we have the perfect solution.

Our QA4 Accessory Stand is a compact, mobile storage system which helps to keep a busy theatre environment tidy and uncluttered.

| Description | Code |
|---|-------|
| QA4 Accessory Stand – W730 x D720 x H1375mm | 21370 |





Gel technology to minimise the risk of pressure sores and infection


Pressure, the amount of force exerted on a given area is often measured in millimetres of mercury (mmHg).

When a force or pressure greater than normal capillary pressure is exerted on a body over time, this can restrict blood flow to the area and cause serious tissue damage.

It is therefore important to reduce this pressure, particularly with the more susceptible parts of the body such as the head in periods of immobility, such as during surgery. Our gel pads help relieve pressure by increasing the area of contact between the body and the supporting surface.


- Latex free
- Reusable






Head Protector


| Description | Code |
|---|-------|
| Azure Head Protector (AZ100) - Infant / Adolescent / Adult - 260 x 250 x 50mm | 10885 |






Prone Head Rest


| Description | Code |
|--|---------|
| Oasis Elite Prone Head Rest (EL141) - Adult - 280 x 240 x 140mm | 10871 |
| Oasis Elite Prone Head Rest (EL142) - Adolescent - 273 x 220 x 135mm | ✂ 10872 |





Supine Head Rest

| Description | Code |
|--|---------|
| Oasis Elite Supine Head Rest (EL151) - Adult - 275 x 237 x 90mm | 10873 |
| Oasis Elite Supine Head Rest (EL152) - Adolescent - 275 x 230 x 65mm | ✂ 10874 |






Face Pad

| Description | Code |
|---|-------|
| Oasis Face Pad (OA140) - 165 x 80 x 8mm | 10849 |









Head and Neck Support


| Description | Code |
|---|-------|
| Oasis+ Head & Neck Support (OP151) - Adult - 264 x 220 x 64mm | 10900 |






Oasis Closed Head Rings


| Description | Code |
|---|-------|
| Oasis Closed Head Ring (OA010) - Neonate / Newborn - Ø90 x 20mm | 10801 |
| Oasis Closed Headring (OA011) - Infant / Adolescent - Ø140 x 35mm | 10802 |
| Oasis Closed Headring (OA012) - Adolescent / Adult - Ø200 x 45mm | 10803 |






Elite Closed Head Rings


| Description | Code |
|--|-------|
| Oasis Elite Closed Head Ring (EL011) - Infant / Adolescent - Ø140 x 34mm | 10861 |
| Oasis Elite Closed Head Ring (EL012) - Adult - Ø200 x 50mm | 10862 |





Open Head Rings

| Description | Code |
|--|-------|
| Oasis Elite Open Head Ring (EL021) - Infant / Adolescent - Ø140 x 34mm | 10863 |
| Oasis Elite Open Head Ring (EL022) - Adult - Ø200 x 50mm | 10864 |





ARM AND HAND

Armboards and Protectors

Flat armboard pads have a slight indentation to cradle the arm and provide effective pressure relief during various surgical procedures.




- Medical grade silicone
- X-ray translucent





Ulnar / Brachial Nerve Protector

| Description | Code |
|---|-------|
| Oasis Ulnar / Brachial Nerve Protector (OA061) - Adolescent / Adult - 460 x 210 x 8mm | 10829 |



Arm Retainer Pad

| Description | Code |
|---|-------|
| Oasis Arm Retainer Pad (OA100) - 170 x 110 x 18mm | 10834 |



Lateral Park Bench Positioner

| Description | Code |
|---|-------|
| Oasis Elite Lateral Park Bench Positioner (EL913) | 10995 |



Oasis Extremity Belt

| Description | Code |
|---|-------|
| Oasis Extremity Belt (OA113) - Pair - 315 x 45 x 10mm | 10843 |




ARM AND HAND



Arm Protector

| Description | Code |
|--|-------|
| Oasis Elite Contoured Armboard / Knee Pad (EL401) - Adult - 510 x 155 x 40mm | 10869 |



Armboard Pads

| Description | Code |
|---|-------|
| Oasis Armboard Pad (OA070) - Infant / Adolescent / Adult - 395 x 115 x 10mm | 10830 |
| Oasis Armboard Pad (OA071) - Infant / Adolescent / Adult - 600 x 200 x 10mm | 10831 |
| Oasis Armboard Pad (OA072) - Adolescent / Adult - 600 x 125 x 10mm | 10832 |



Arm Protector

| Description | Code |
|---|-------|
| Azure Arm Protector (AZ201) - Adolescent / Adult - 460 x 180 x 35mm | 10890 |






Protection for the wrist and elbow, particularly the ulnar nerve

LEG, KNEE AND FOOT

Legs and Heels


The anatomy of the heel – its combination of nerves and lack of soft tissue over the bone – make it one of the most vulnerable parts of the body during surgery, particularly when the patient is in a supine position.

- Easily Cleaned
- High load bearing capacity



Heel Cup

| Description | Code |
|---|-------|
| Oasis Heel Cup (OA131) - Neonate / Infant / Adolescent / Adult - Pair - 105 x 90 x 43mm | 10846 |




Contoured Knee / Arm Pads

| Description | Code |
|---|-------|
| Oasis Contoured Armboard / Knee Pad (OA101) - 510 x 150 x 35mm | 10833 |
| Oasis Elite Contoured Armboard / Knee Pad (EL402) - Adolescent / Adult - 395 x 140 x 40mm | 10870 |



Dome Positioner

| Description | Code |
|---|-------|
| Oasis Elite Dome Positioner (EL210) - Adolescent / Adult - 350 x 150 x 75mm | 10876 |



Lateral Leg Positioner

| Description | Code |
|--|-------|
| Oasis Elite Lateral Leg Positioner (EL800) - Adult - 650 x 400 x 230mm | 10992 |



LEG, KNEE AND FOOT



Heel Pad

| Description | Code |
|---|-------|
| Oasis Elite Heel Pad (EL130) - Adolescent / Adult - Pair - 180 x 100 x 70mm | 10875 |



Heel Protector

| Description | Code |
|---|-------|
| Azure Heel Protector (AZ401) - Infant / Adolescent / Adult - 190 x 175 x 30mm | 10891 |



Leg Protector

| Description | Code |
|---|-------|
| Azure Leg Protector (AZ400) - Adolescent / Adult - 515 x 180 x 40mm | 10892 |



Heel and Ankle Protector

| Description | Code |
|--|-------|
| Oasis Heel & Ankle Protector (OA200) - Neonate / Infant / Adolescent / Adult - 1250 x 120 x 10mm | 10853 |



Protective pressure care for vulnerable extremities like the heel



DESIGN & INNOVATION



TECHNICAL EXCELLENCE



GLOBAL DISTRIBUTION

SHOULDER, SPINE AND TORSO


Trulife Oasis Elite option:
light weight and pre-contoured

One of four pressure care construction options from Trulife, the Oasis Elite range combines silicone gel and light weight foam with a pre-contoured surface to provide superior pressure relief and patient protection in theatre.

Anatomical pre-contouring facilitates greater conformity between gel and body shape; weight is then redistributed over a larger area.

Oasis Elite products are more than 50% lighter than comparable gel-only products.

- Foam and silicone gel combination
- Enhanced protection
- Light weight
- Highly durable
- Easy to clean and maintain



Montreal Laminectomy Support

| Description | Code |
|---|---------|
| Elite Foam Large Prone Positioner (EL904) - Montreal Laminectomy - Adult - 1115 x 520 x 140mm | ✂ 10990 |




Prone Chest Positioner

| Description | Code |
|--|---------|
| Oasis Elite Prone Positioner (EL903) - Adult - 610 x 530 x 130mm | ✂ 10883 |



Prone Positioner

| Description | Code |
|--------------------------------------|---------|
| Oasis Elite Prone Positioner (EL905) | ✂ 10884 |



Elite Prone Positioning System

| Description | Code |
|---|------------|
| Elite Foam ARDS Prone Positioning System (ELP200) | ✂ 10991 |
| Elite Foam Disposable Prone Head Positioner (ELP141) Supplied in packs of 3 | ✂ 10991-A1 |



SHOULDER, SPINE AND TORSO




Oasis Chest Rolls

| Description | Code |
|--|-------|
| Oasis Flat Bottom Chest Roll (OA211) - Adult - 305 x 95 x 75mm | 10854 |
| Oasis Flat Bottom Chest Roll (OA212) - Adult - 405 x 95 x 75mm | 10855 |
| Oasis Flat Bottom Chest Roll (OA213) - Adult - 505 x 92 x 75mm | 10856 |



Elite Chest Rolls

| Description | Code |
|--|---------|
| Oasis Elite Flat Bottom Chest Roll (EL214) - Adult - 300 x 140 x 110mm | ✂ 10877 |
| Oasis Elite Flat Bottom Chest Roll (EL215) - Adult - 400 x 140 x 110mm | ✂ 10878 |
| Oasis Elite Flat Bottom Chest Roll (EL216) - Adult - 505 x 140 x 110mm | ✂ 10879 |
| Oasis Elite Flat Bottom Chest Roll (EL217) - Adult - 405 x 150 x 150mm | ✂ 10880 |
| Oasis Elite Flat Bottom Chest Roll (EL218) - Adult - 510 x 150 x 150mm | 10881 |



Patient Positioner

| Description | Code |
|--|---------|
| Azure Universal Positioner (AZ500) - Adolescent / Adult - 425 x 155 x 60mm | ✂ 10893 |



DESIGN & INNOVATION



TECHNICAL EXCELLENCE



GLOBAL DISTRIBUTION

Secure supine, lateral or prone positioning

Trulife Azure option:
redefining pressure care


Azure’s revolutionary gel combination works together to replicate and act as an additional layer of protective tissue and skin.

It conforms to the individual's shape in a more contoured and less structured manner. This results in effective protection against nerve damage as well as pressure sore prevention.

Azure: advanced construction for state-of-the-art pressure care




- Advanced surgical positioning and protection
- Revolutionary new gel combination
- Contoured protection
- Redefining pressure care
- 12 month guarantee



Crutch Stirrup Pads

| Description | Code |
|---|-------|
| Oasis Crutch Stirrup Pad (OA090) - Adolescent / Adult - Pair - 320 x 220 x 10mm | 10838 |
| Oasis Crutch Stirrup Pad (OA091) - Adolescent / Adult - Pair - 395 x 245 x 10mm | 10839 |



Crutch Stirrup Pad

| Description | Code |
|---|-------|
| Oasis Crutch Stirrup Pad (OA092) - Adolescent / Adult - Pair - 520 x 300 x 10mm | 10840 |




Table Pads

| Description | Code |
|---|---------|
| Oasis 3-piece Operating Table Pad (OA037) - Torso Section with Cut-out - 900 x 500 x 10mm | 10813 |
| Oasis Perineal Table Pad (OA040) - 520 x 520 x 10mm | ✂ 10816 |



Profile Pad – Universal

| Description | Code |
|--|---------|
| Azure Profile Pad (AZ610) - Universal - 520 x 520 x 10mm | ✂ 10886 |

Profile Pad – ¾ Length

| Description | Code |
|--|---------|
| Azure Profile Pad (AZ611) - ¾ Length - 1150 x 520 x 10mm | ✂ 10887 |



Sacral Protector

| Description | Code |
|--|-------|
| Azure Sacral Protector (AZ300) - Adolescent / Adult - 400 x 380 x 25mm | 10888 |



Versatile pressure care pads for a variety of procedures



DESIGN & INNOVATION



TECHNICAL EXCELLENCE



GLOBAL DISTRIBUTION

Trulife Oasis option:
the original silicone solution

One of four pressure care construction options from Trulife, the Oasis range consists of reusable 100% silicone gel pads.

The silicone is softer than both skin and underlying tissue which means that pressure is relieved by the slight movement of the silicone gel and resulting dissipation of pressure/force across the product.

Oasis: gel pressure care for positioning across a wide range of surgical specialisms



- Silicone superior pressure relief
- Conforms to body shape
- Comprehensive range
- Can be heated or cooled
- Latex free and skin compatible





Patient Positioning

| Description | Code |
|--|-------|
| Oasis Universal Patient Positioner (OA050) - Infant / Adolescent / Adult - 150 x 60 x 40mm | 10823 |
| Oasis Universal Patient Positioner (OA051) - Infant / Adolescent / Adult - 330 x 70 x 50mm | 10824 |
| Oasis Universal Patient Positioner (OA052) - Adolescent / Adult - 290 x 100 x 60mm | 10825 |
| Oasis Universal Patient Positioner (OA053) - Adolescent / Adult - 345 x 120 x 70mm | 10826 |



Operating Table Pads

| Description | Code |
|---|-------|
| Oasis Standard Operating Table Pad (OA030) - Universal - 520 x 520 x 10mm | 10807 |




Operating Table Pads

| Description | Code |
|---|-------|
| Full Body - Oasis Standard Operating Table Pad (OA031) - ¾ Length - 1150 x 520 x 10mm | 10808 |
| Oasis Standard Operating Table Pad (OA032) - Full Length - 1800 x 520 x 10mm | 10809 |



Oasis Patient Positioning Belt

| Description | Code |
|---|-------|
| Oasis Patient Positioning Belt (OA110) - 505 x 85 x 8mm | 10842 |



Oasis Non-slip Discs

| Description | Code |
|---|-------|
| Oasis+ Non-slip Disc (OP015) - Ø200 x 4mm | 10975 |



Vacuum Beanbag Positioners

| Description | Code |
|---|-------|
| Vacuum Beanbag Positioner - Full Body - 1000 x 900mm | 10275 |
| Vacuum Beanbag Positioner - Square - 500 x 500mm | 10276 |
| Vacuum Beanbag Positioner - Head Stabiliser - 400 x 300mm | 10278 |
| Vacuum Beanbag Positioner - Full Length - 1650 x 720mm | 10279 |



PATIENT POSITIONS

PURCHASE GUIDE

SUPINE position

- Administration of General Anaesthesia
- Abdominal surgery
- Open heart surgery
- Surgery of the face, neck or mouth
- Most surgery on extremities



The patient is on their back with their face up towards the ceiling.
Their legs are uncrossed and their arms are at their sides, or on arm boards.

(Includes Trendelenburg and reverse Trendelenburg)

PRONE position

- Surgeries to the back, spine and the backs of the legs



The patient is on their abdomen and their arms are at their sides or on armboards.

(Includes Kraske 'Jack Knife' position)

LATERAL position

- Surgery of the chest or lungs
- Surgery on the kidneys
- Hip surgery



The patient is on their side, with their bottom leg flexed and their top leg straight.

The bottom arm is supported on an arm board, and the top arm has special support.

LITHOTOMY position

- Gynaecological surgery
- Genitourinary (GU) surgery



The patient is on their back, their legs are lifted in stirrups and their arms are usually on arm boards.

| Code | Description | POSITIONS | | | |
|-------|--|-----------|-------|---------|-----------|
| | | SUPINE | PRONE | LATERAL | LITHOTOMY |
| 10802 | Oasis Closed Headring (OA011) - Infant / Adolescent - Ø140 x 35mm | ✓ | | | ✓ |
| 10803 | Oasis Closed Headring (OA012) - Adolescent / Adult - Ø200 x 45mm | ✓ | | | ✓ |
| 10807 | Oasis Standard Operating Table Pad (OA030) - Universal - 520 x 520 x 10mm | ✓ | ✓ | ✓ | ✓ |
| 10808 | Oasis Standard Operating Table Pad (OA031) - ¾ Length - 1150 x 520 x 10mm | ✓ | ✓ | ✓ | ✓ |
| 10809 | Oasis Standard Operating Table Pad (OA032) - Full Length - 1800 x 520 x 10mm | ✓ | | | ✓ |
| 10813 | Oasis 3-piece Operating Table Pad (OA037) - Torso Section with Cut-out - 900 x 500 x 10mm | ✓ | | | ✓ |
| 10824 | Oasis Universal Patient Positioner (OA051) - Infant / Adolescent / Adult - 330 x 70 x 50mm | ✓ | ✓ | ✓ | |
| 10825 | Oasis Universal Patient Positioner (OA052) - Adolescent / Adult - 290 x 100 x 60mm | ✓ | ✓ | ✓ | |
| 10826 | Oasis Universal Patient Positioner (OA053) - Adolescent / Adult - 345 x 120 x 70mm | ✓ | ✓ | ✓ | |
| 10829 | Oasis Ulnar / Brachial Nerve Protector (OA061) - Adolescent / Adult - 460 x 210 x 8mm | ✓ | ✓ | ✓ | ✓ |
| 10831 | Oasis Armboard Pad (OA071) - Infant / Adolescent / Adult - 600 x 200 x 10mm | ✓ | ✓ | ✓ | ✓ |
| 10838 | Oasis Crutch Stirrup Pad (OA090) - Adolescent / Adult - Pair - 320 x 220 x 10mm | ✓ | ✓ | ✓ | ✓ |
| 10839 | Oasis Crutch Stirrup Pad (OA091) - Adolescent / Adult - Pair - 395 x 245 x 10mm | ✓ | ✓ | ✓ | ✓ |
| 10853 | Oasis Heel & Ankle Protector (OA200) - Neonate / Infant / Adolescent / Adult - 1250 x 120 x 10mm | ✓ | | ✓ | |
| 10856 | Oasis Flat Bottom Chest Roll (OA213) - Adult - 505 x 92 x 75mm | ✓ | ✓ | ✓ | |
| 10861 | Oasis Elite Closed Head Ring (EL011) - Infant / Adolescent - Ø140 x 34mm | ✓ | | | ✓ |
| 10862 | Oasis Elite Closed Head Ring (EL012) - Adult - Ø200 x 50mm | ✓ | | | ✓ |
| 10863 | Oasis Elite Open Head Ring (EL021) - Infant / Adolescent - Ø140 x 34mm | | ✓ | ✓ | |
| 10864 | Oasis Elite Open Head Ring (EL022) - Adult - Ø200 x 50mm | | ✓ | ✓ | |
| 10869 | Oasis Elite Contoured Armboard / Knee Pad (EL401) - Adult - 510 x 155 x 40mm | ✓ | ✓ | ✓ | ✓ |
| 10871 | Oasis Elite Prone Head Rest (EL141) - Adult - 280 x 240 x 140mm | | ✓ | | |
| 10873 | Oasis Elite Supine Head Rest (EL151) - Adult - 275 x 237 x 90mm | ✓ | | | |
| 10875 | Oasis Elite Heel Pad (EL130) - Adolescent / Adult - Pair - 180 x 100 x 70mm | ✓ | | | |
| 10876 | Oasis Elite Dome Positioner (EL210) - Adolescent / Adult - 350 x 150 x 75mm | | ✓ | ✓ | |
| 10881 | Oasis Elite Flat Bottom Chest Roll (EL218) - Adult - 510 x 150 x 150mm | ✓ | ✓ | ✓ | |
| 10882 | Oasis Elite Lateral Positioner (EL900) - Adolescent / Adult - 720 x 500 x 160mm | | ✓ | | |
| 10884 | Oasis Elite Prone Positioner (EL905) - 675 x 500 x 100mm | | ✓ | | |
| 10885 | Azure Head Protector (AZ100) - Infant / Adolescent / Adult - 260 x 250 x 50mm | ✓ | | | ✓ |
| 10886 | Azure Profile Pad (AZ610) - Universal - 520 x 520 x 10mm | ✓ | ✓ | ✓ | ✓ |
| 10887 | Azure Profile Pad (AZ611) - ¾ Length - 1150 x 520 x 10mm | ✓ | ✓ | ✓ | ✓ |
| 10888 | Azure Sacral Protector (AZ300) - Adolescent / Adult - 400 x 380 x 25mm | | | | ✓ |
| 10890 | Azure Arm Protector (AZ201) - Adolescent / Adult - 460 x 180 x 35mm | ✓ | ✓ | ✓ | ✓ |
| 10891 | Azure Heel Protector (AZ401) - Infant / Adolescent / Adult - 190 x 175 x 30mm | ✓ | | | |
| 10892 | Azure Leg Protector (AZ400) - Adolescent / Adult - 515 x 180 x 40mm | ✓ | | | |
| 10990 | Elite Foam Large Prone Positioner (EL904) - Montreal Laminectomy - Adult - 1115 x 520 x 140mm | | ✓ | | |
| 10992 | Oasis Elite Lateral Leg Positioner (EL800) - Adult - 650 x 400 x 230mm | | | ✓ | |



Scan for Trulife pressure care training video



DESIGN & INNOVATION



TECHNICAL EXCELLENCE



GLOBAL DISTRIBUTION

INSTRUMENT TROLLEYS

We offer a wide range of Instrument Trolleys and every component is quality hallmarked.

Each product is individually inspected and available for immediate dispatch.

All our Instrument Trolleys have Fixed Downturned Shelves, come complete with Stainless Steel bright polished frames and use only British Standard recommended gauge 304 stainless steel. All shelves are 2.0mm thick with a 240 silicon grain ultra smooth satin finish. The castors are twin-wheeled, non-corrosive and ideal for modern clinical areas.



| Description | Code |
|---|-------|
| Instrument Trolley W450 x D450 x H850mm – Fixed Downturned Shelves | 25000 |
| Instrument Trolley W600 x D450 x H850mm – Fixed Downturned Shelves | 25010 |
| Instrument Trolley W750 x D450 x H850mm – Fixed Downturned Shelves | 25020 |
| Instrument Trolley W900 x D450 x H850mm – Fixed Downturned Shelves | 25030 |
| Instrument Trolley W1200 x D450 x H850mm – Fixed Downturned Shelves | 25040 |
| Instrument Trolley W1200 x D600 x H850mm – Fixed Downturned Shelves | 25050 |



Does your equipment migrate to other departments?

Use our new Coloured Castors to easily identify where individual pieces of equipment belong.

Available as a factory fit in Red, Blue, Green, Yellow, Purple, Black or White.



- High strength plastic-bodied castors do not suffer corrosion like traditional steel-bodied models and are quick and simple to replace if damaged

DRESSING / EQUIPMENT TROLLEYS



Dressing trolleys have Fixed Upturned Shelves which are ideal for anaesthetic and recovery room areas.

Optional swivel bowls for plaster rooms are available on request.

General Specifications: Instrument and Dressing / Equipment Trolleys

| | |
|-----------------|--|
| FRAMES | Stainless Steel bright polished tube |
| SHELVES | Stainless Steel polished sheet 240 silicon grain ultra smooth satin finish |
| CASTORS | 65mm Twin-wheeled, polyurethane-tyred, non-corrosive |
| WEIGHT CAPACITY | 100kg |

| Description | Code |
|--|-------|
| Dressing Trolley W450 x D450 x H850mm – Fixed Upturned Shelves | 25100 |
| Dressing Trolley W750 x D450 x H850mm – Fixed Upturned Shelves | 25120 |
| Dressing Trolley W900 x D450 x H850mm – Fixed Upturned Shelves | 25130 |

Optional Features

| Description | Code |
|--|-------|
| Theatre Furniture – Castors – Ø65mm Braking Rotating Buffers | 25367 |
| Theatre Furniture – Castors – Red | 25290 |
| Theatre Furniture – Castors – Blue | 25220 |
| Theatre Furniture – Castors – Green | 25222 |
| Theatre Furniture – Castors – Yellow | 25224 |
| Theatre Furniture – Castors – Purple | 25226 |
| Theatre Furniture – Castors – Black | 25227 |
| Theatre Furniture – Castors – White | 25228 |
| Factory fit – please specify at time of ordering | 25229 |

STAINLESS STEEL FURNITURE



DESIGN & INNOVATION



TECHNICAL EXCELLENCE



GLOBAL DISTRIBUTION

SURGEON'S SEATING

Our gas-lift stools are made in a modular design which allows customers to choose the exact specifications to meet their needs.



General Specifications: Surgeon's Stools

| | |
|-----------------|---|
| BASE | 5-Star aluminium Teflon® coated finish which is ultra clean and wear resistant – 585mm diameter |
| TOP | Polished Stainless Steel |
| HEIGHT RANGE | 480mm – 680mm |
| CASTORS | 65mm Twin-wheeled, polyurethane-tyred, non-corrosive |
| SEAT | 380mm Pressure care pad |
| BACKREST | Adjustable, polished Stainless Steel with pad |
| ARMRESTS | Polished Stainless Steel with pad |
| WEIGHT CAPACITY | 150kg |

Optional Features

| Description | Code |
|--|-------|
| Surgeon's Stool – Backrest – includes Pad | 25320 |
| Surgeon's Stool – Armrests – includes Pads | 25330 |
| Surgeon's Stool – Castors – Ø65mm Twin Wheel – Braking | 25370 |
| Surgeon's Stool – Castors – Ø65mm Twin Wheel – Brake-when-loaded | 25375 |
| Theatre Furniture – Castors – Red | 25220 |
| Theatre Furniture – Castors – Blue | 25222 |
| Theatre Furniture – Castors – Green | 25224 |
| Theatre Furniture – Castors – Yellow | 25226 |
| Theatre Furniture – Castors – Purple | 25227 |
| Theatre Furniture – Castors – Black | 25228 |
| Theatre Furniture – Castors – White | 25229 |
| Factory fit – please specify at time of ordering | |



- Stool mechanism enclosed within moulded cover for easy cleaning

Transfusion / Infusion Stand 5-Star Base

Transfusion / Infusion Stand Weighted 5-Star Base

Transfusion / Infusion Stand Ward Environment



Stools are foot-operated, which allows adjustment without compromising the sterile field. We have maintained a simple form by enclosing the mechanism inside a moulded cover and all metalwork is mirror polished with a specialist heavy duty Teflon® coated finish on the base. Pressure Care pads are fitted as standard.

- Foot-operated surgeon's stools allow adjustment without compromising the sterile field



- Modular design enables customers to tailor each stool to their exact requirements

| Description | Code |
|--|-------|
| Surgeon's Stool – Castors – includes Seat Cushion | 25300 |
| Surgeon's Stool – Stump Feet – includes Seat Cushion | 25310 |

TRANSFUSION / INFUSION STANDS – MOBILE

We supply three types of Transfusion Stand.

We offer a light weight, 5-Star Based version with superb mobility, a Heavy Cast Iron Based version with exceptional stability, PLUS – a special budget version is now available for ward environments.

- Provita easy action single handed mechanism: quick release, safety locking†

- Extendable and removable, unique hook mechanism



General Specifications: Transfusion / Infusion Poles

| | | |
|------------------|-----------------|---|
| 5-Star Base | BASE | 5-Star aluminium Teflon® coated finish which is ultra clean and wear resistant – 585mm diameter |
| | POST / COLUMN | Bright polished Stainless Steel |
| | CASTORS | 65mm Polyurethane-tyred, non-corrosive |
| | HEIGHT RANGE | 1505mm – 2185mm |
| | WEIGHT | 5kg |
| Heavy Base | WEIGHT CAPACITY | 12kg |
| | BASE | Resin filled acrylic capped ABS vac-formed cover, antibacterial – 570mm diameter |
| | POST / COLUMN | Bright polished Stainless Steel |
| | CASTORS | 65mm Polyurethane-tyred, non-corrosive |
| | HEIGHT RANGE | 1540mm – 2220mm |
| Ward Environment | WEIGHT | 17kg |
| | WEIGHT CAPACITY | 25kg |
| | BASE | Polypropylene |
| | POST / COLUMN | Bright polished Stainless Steel |
| | CASTORS | 65mm Polyurethane-tyred, non-corrosive |
| | HEIGHT RANGE | 1285mm – 2110mm |
| | WEIGHT | 3.5kg |
| | WEIGHT CAPACITY | 4kg |

Our transfusion stands incorporate several important features rarely found on other transfusion stands on the market. These include solid stainless steel poles, an innovative quick release safety locking mechanism† and drainage at the base of the column to prevent fluid build-up. There is a specialist heavy duty Teflon® coated finish on the light weight base† and a unique acrylic capped antibacterial cover to the heavy base version.

| Description | Code |
|---|-------|
| Transfusion / Infusion Stand – Mobile – 5-Star Base – 4 Hook | 25420 |
| Transfusion / Infusion Stand – Mobile – Weighted 5-Star Base – 4 Hook | 25430 |
| Transfusion / Infusion Stand – Mobile – 5-Star Base – 4 Hook (Ward Environment) | 25440 |

Optional Features

| Description | Code |
|--|-------|
| Theatre Furniture – Castors – Red | 25220 |
| Theatre Furniture – Castors – Blue | 25222 |
| Theatre Furniture – Castors – Green | 25224 |
| Theatre Furniture – Castors – Yellow | 25226 |
| Theatre Furniture – Castors – Purple | 25227 |
| Theatre Furniture – Castors – Black | 25228 |
| Theatre Furniture – Castors – White | 25229 |
| Factory fit – please specify at time of ordering | |

†Excludes Ward Environment version



Designing and manufacturing products with infection control built in

MAYO TABLES AND BOWL STANDS

Mayo Tables

Supplied with two base options, both weighted for exceptional stability, Mayo tables are manufactured to a traditional design and recognised for their ease of use and functionality.



General Specifications: Mayo Tables

| | |
|-----------------------|--|
| FRAME | Stainless Steel electropolished finish |
| CASTORS | 65mm Twin-wheeled, polyurethane-tyred, non-corrosive |
| HEIGHT RANGE | 870mm – 1220mm |
| BASE WIDTH | 495mm |
| BASE LENGTH | 670mm |
| HEIGHT TO TOP OF BASE | 225mm |
| TRAY SIZE | 510 x 380 x 25mm |

| Description | Code |
|---|-------|
| Mayo Table – 4 Castors | 25500 |
| Mayo Table – 2 Castors & 2 Flattened Forefeet | 25510 |

Bowl Stands

Manoeuvrable and functional, these models are manufactured with a Teflon® coated 5-Star Base to ensure maximum stability.



Bowl/s included

Optional Features

| Description | Code |
|--|-------|
| Theatre Furniture – Castors – Red | 25220 |
| Theatre Furniture – Castors – Blue | 25222 |
| Theatre Furniture – Castors – Green | 25224 |
| Theatre Furniture – Castors – Yellow | 25226 |
| Theatre Furniture – Castors – Purple | 25227 |
| Theatre Furniture – Castors – Black | 25228 |
| Theatre Furniture – Castors – White | 25229 |
| ☞ Factory fit – please specify at time of ordering | |

General Specifications: Bowl Stands

| | |
|------------|---|
| FRAME | Stainless Steel electropolished finish |
| BASE | 5-Star aluminium Teflon® coated finish which is ultra clean and wear resistant – 585mm diameter |
| BOWLS | 350mm Stainless Steel electropolished |
| CASTORS | 65mm Twin-wheeled, polyurethane-tyred, non-corrosive |
| BUFFERING | Spiral buffers on rings to eliminate rattle |
| DIMENSIONS | Height / Base Diameter / Width |
| Single | 860mm / 585mm / 585mm |
| Double | 860mm / 585mm / 755mm |

| Description | Code |
|---|-------|
| Bowl Stand – Single Bowl | 25600 |
| Bowl Stand – Double Bowl – Side-by-Side Mounted | 25610 |

KICK-ABOUT AND SURGEON'S PLATFORM

Kick-about

These models are made with the highest grade materials. All are electropolished to produce an even smoother surface and finished with twin-wheeled, polyurethane-tyred, non-corrosive castors.




General Specifications: Kick-about

| | |
|-----------------|--|
| FRAME | Stainless Steel electropolished finish |
| BOWL | Stainless Steel electropolished and graduated finish |
| CASTORS | 65mm Twin-wheeled, polyurethane-tyred, non-corrosive |
| BUFFERING | Spiral buffers on rings to eliminate rattle |
| DIMENSIONS | Height / Diameter |
| Kick-about Bowl | 230mm / 370mm |

| Description | Code |
|-------------------------------------|-------|
| Kick-about Bowl | 25700 |
| Kick-about Bucket – 12 Litre Bucket | 25710 |

Surgeon's Platform

Designed to allow full surgical access, our platforms achieve maximum stability by using only 4.5mm non-slip aluminium treadplate.



General Specifications: Surgeon's Platform

| | |
|------------|--|
| FRAME | Stainless Steel electropolished finish |
| PLATFORM | Aluminium 4.5mm non-slip treadplate |
| DIMENSIONS | Height / Width / Length |
| Frame | 190mm / 300mm / 760mm |
| Treadplate | 190mm / 300mm / 610mm |

| Description | Code |
|--------------------|-------|
| Surgeon's Platform | 25800 |

Linen Receiver

This stable but highly manoeuvrable item has a mild steel powder coated frame and twin-wheeled, polyurethane-tyred, non-corrosive castors.



General Specifications: Linen Receiver

| | |
|------------|--|
| FRAME | Mild steel powder coated in white |
| CASTORS | 65mm Twin-wheeled, polyurethane-tyred, non-corrosive |
| BAG | White terylene includes draw cord |
| DIMENSIONS | Height / Ring Diameter / Width |
| Single | 870mm / 490mm / 490mm |

| Description | Code |
|---|-------|
| Soiled Linen Receiver – Single Receiver – includes Terylene Bag | 25900 |

Practical products: manoeuvrable, stable and functional



DESIGN & INNOVATION



TECHNICAL EXCELLENCE



GLOBAL DISTRIBUTION

GENERAL

Spare components



| Description | Code |
|--|--------|
| Surgeon's Stool – Cushion – Ø380mm | 25340 |
| Surgeon's Stool – Backrest Pad | 25350 |
| Surgeon's Stool – Armrest Pad | 25360 |
| Mayo Table – Stainless Steel Tray Only | 25515 |
| Bowl – Stainless Steel – Ø350mm – H140mm – 6 Litre Capacity | 25705 |
| Castor – Ø65mm Twin Wheel – 5-star Base – Grey – Braking | 692005 |
| Castor – Ø65mm Twin Wheel – Surgeon's Stool – Grey – Brake-when-loaded | 692010 |
| Castor – Ø65mm Twin Wheel – Tubular Leg – Grey | 692015 |
| Castor – Ø65mm Twin Wheel – Tubular Leg – Grey – Braking | 692020 |
| Castor – Ø65mm Twin Wheel – 5-star Base & Tubular Leg – Red | 692025 |

| Description | Code |
|--|--------|
| Castor – Ø65mm Twin Wheel – 5-star Base & Tubular Leg – Blue | 692030 |
| Castor – Ø65mm Twin Wheel – 5-star Base & Tubular Leg – Green | 692035 |
| Castor – Ø65mm Twin Wheel – 5-star Base & Tubular Leg – Yellow | 692040 |
| Castor – Ø65mm Twin Wheel – 5-star Base & Tubular Leg – Purple | 692045 |
| Castor – Ø65mm Twin Wheel – 5-star Base & Tubular Leg – Black | 692050 |
| Castor – Ø65mm Twin Wheel – 5-star Base & Tubular Leg – White | 692055 |
| Castor Insert – for ID Ø23mm Tubular Leg | 692060 |
| Surgeon's Stool – Stump Foot | 692075 |
| Rubber Ferrules / Feet – Ø19mm Tubular Leg | 692080 |
| Rubber Ferrules / Feet – Ø25mm Tubular Leg | 692085 |

Touch Fastening



| Description | Code |
|--|--------|
| Touch Fastener Pack - includes 2 x 500mm Hook & Loop Strips plus 12 x VELCRO® Brand Ø45mm Loop Coins | 691000 |
| VELCRO® Brand Touch Fastener - Ø45mm Loop Coin | 691005 |

Rubber Ferrules



Available in two different sizes, we offer non-slip replacement rubber ferrules for items such as surgeon's platforms.

| Description | Code |
|--|--------|
| Theatre Furniture - Rubber Ferrules / Feet - for Ø19mm Tubular Leg | 692080 |
| Theatre Furniture - Rubber Ferrules / Feet - for Ø25mm Tubular Leg | 692085 |

Storage Cart for Lift-Assist™ Stirrups

The perfect way to store and protect your investment, this strong and light weight cart holds and transports stirrups or leg holders.

- Keeping busy theatre areas tidy and uncluttered



General Specifications: Lift-Assist Stirrups Stand

| | |
|-----------------------|--|
| FRAME | Roto-moulded liner, low-density polyethylene |
| HANGERS / BASE COVERS | Roto-moulded liner, low-density polyethylene |
| CASTORS | 50mm Twin-wheeled, nylon |
| DIMENSIONS | Height / Width / Depth |
| | 1370mm / 530mm / 530mm |

| Description | Code |
|--|-------|
| Storage Cart – for Lift-Assist™ Lithotomy Stirrups | 10760 |

QA3 & QA4 STRETCHERS

QA3 Patient, Ophthalmic and Emergency Stretchers

Ideal for patient transport, examination and intubation in theatre and clinical areas – with specially designed versions for Ophthalmics, Accident and Emergency.



| Description | Code |
|---|-------|
| QA3™ Patient Stretcher – includes K8 Pressure Care Mattress | 21110 |
| QA3 DRIVE™ Patient Stretcher – includes K8 Pressure Care Mattress | 21112 |
| QA3 DRIVE™ Powered Patient Stretcher – includes K8 Pressure Care Mattress | 21114 |
| QA3™ Powered Patient Stretcher – includes K8 Pressure Care Mattress | 21116 |
| QA3™ Ophthalmic Stretcher – includes K8 Pressure Care Mattress | 21118 |
| QA3™ Emergency Stretcher – with Radiolucent Platform – includes K8 Pressure Care Mattress | 21120 |
| QA3 DRIVE™ Emergency Stretcher – with Radiolucent Platform – includes K8 Pressure Care Mattress | 21122 |
| QA3 DRIVE™ Powered Emergency Stretcher – with Radiolucent Platform – includes K8 Pressure Care Mattress | 21124 |
| QA3™ Powered Emergency Stretcher – with Radiolucent Platform – includes K8 Pressure Care Mattress | 21126 |



QA4 Mobile Surgery System



Whether for Day / Ambulatory Surgery or Mainstream Theatre, the QA4 is the complete multi-purpose mobile operating system.

With five remotely-activated powered functions – raise and lower, back rest, traversing top, lateral tilt and Trendelenburg – the stretcher offers full 'C' arm coverage and superior surgical access. It also offers exceptional stability with a patient weight capacity of 250kg.

| Description | Code |
|--|-------|
| QA4™ Mobile Surgery System | 21300 |
| QA4™ Accessory Stand – W730 x D720 x H1375mm | 21370 |



Stretcher solutions for transport treatment and recovery